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en INSTRUCTIONS FOR USE

Intended use

Used to convey the ventilation gases (air and oxygen) from the source to the interface device during non-invasive ventilation.

Restrictions of use

- The device must be used by qualified and trained medical/nursing staff.
- To be used for administration of air and oxygen.
- If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.

Instructions for use

Remove the device from the package.

Attach the terminal connectors of the tubes to the gas distributor (ventilator/ flow meter/Venturi meter) and to the device on the patient side (hood/mask/ etc.). Check that the connection is secure before proceeding with administration of the therapy gases.

Warnings and precautions

⚠ Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

- In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
- Monitor the clinical parameters of the patient during the treatment.
- It is unadvisable to use humidified ventilation gas since this would cause the formation of condensation on the walls of the device.
- The device can be used for maximum 7 days. Depending on the clinical picture of the patient the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
- Non-sterile. Do not sterilize.
- Do not use on another patient.
- Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).
- Dispose of the materials immediately after use conforming to the current laws and regulations.
- Use only with devices that have connections in compliance with the ISO 5356-1 standards.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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