## LEGENDA SIMBOLI/SYMBOL LEGEND

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**en****Intended use**

Used for oral care by healthcare operators in not self-sufficient patients.

**Indications for use**

Care of the oral cavity and mechanical removal of plaque in order to prevent diseases and infections of the oral cavity and then potential VAPs for not self-sufficient patients.

**Contraindications**

Do not use in case of bleeding in the oral cavity.

**Restrictions of use**

The device must be used by qualified and trained medical/nursing staff.

Male terminal connector which can be connected to suction tubes (connected to suction systems considered as medical devices).

Use in combination with ORAL-FIT mouthwash, available in 15 ml vials, sold separately.

**Complications/side effects**

None if used according to the instructions for use.

**Instructions for use**

1. Unscrew/reassemble the toothbrush, if necessary.
2. Connect the toothbrush to suction tubing using the terminal male connector. Check that the connection is secure before proceeding with suction.
3. Close the control port with thumb to activate suction.
4. Brush gums, teeth, tongue, palate and cheeks.
5. Remove dental plaque residues through suction.
6. Wash toothbrush thoroughly (both inside and outside) with plenty of water. Dry with care, possibly using clean gauze.

**Warnings / Precautions**

In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.

If the patient has no teeth, brush the gums, tongue and oral cavity gently.

Act gently on the tongue to prevent patient's vomit.

Non-sterile. Do not sterilize.

Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).

Dispose of the materials immediately after use CONFORMING TO the current laws and regulations Dispose of the device after use in conformity with the legislation in force as regards hospitals.

**Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**