
<b>it</b> TAPPO CALIBRATO	<b>fr</b> BOUCHON CALIBRE	<b>es</b> TAPÓN CALIBRADO	<b>nl</b> GEKALIBREERDE DOP	<b>el</b> ΒΑΘΜΟΝΟΜΗΜΕΝΟ ΠΩΜΑ
<b>en</b> PRESET CAP	<b>de</b> KALIBRIERTER VERSCHL USS	<b>pt</b> TAMPÁ CALIBRADA	<b>sv</b> REGLERTAPP	<b>sl</b> KALIBRIRANI ZAMAŠEK
<b>sk</b> KALIBROVANÝ UZÁVER	<b>tr</b> AYARLI KAPAK	<b>sr</b> KALIBRIRAN POKLOPAC		

Istruzioni per l'uso / Instructions for use / Mode d'emploi / Gebrauchsanleitung / Instrucciones de uso / Instruções de utilização / Gebruiksaanwijzing / Bruksanvisning / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / Navodilo za uporabo / Návod na použitie / Kullanım talimatları / Uputstvo za upotrebu

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**en** INTENDED USE

The PEEP (Positive End Expiratory Pressure) valve is used in the respiratory systems to maintain a positive end expiratory pressure in the airways.

**INSTRUCTIONS FOR USE**

Remove the device from the package and connect it to the outlet connector of the helmet. The device is now ready for operation. Activate gas flow.

**WARNINGS**

Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients

- In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
- The device must be used by qualified medical/nursing staff aware of the possible hemodynamic side effects for the patient caused by positive end respiratory pressure treatments.
- Monitor the clinical parameters of the patient during the treatment and always check the ventilation pressure value with a pressure gauge connected to a pressure port on the circuit.
- The compliance of the device to the above pressure/flow table is approximate
- The device can be used for maximum 7 days. Depending on the clinical condition of the patient the medical/nursing staffs are responsible for an eventual more frequent replacement of the device.
- Non-sterile. Do not sterilize.
- Do not use on another patient.
- Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).
- Dispose of the materials immediately after use CONFORMING TO the current laws and regulations.
- Use only with devices that have connections in compliance with the ISO 5356-1 standards.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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