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en INSTRUCTIONS FOR USE

INTENDED USE

Pressure relief valve for helmet, suitable to prevent overpressure in positive pressure respiratory system and the risk of barotraumas.

RESTRICTIONS OF USE

The device must be used by qualified and trained medical/nursing staff.

If used properly, the device may be used continuously for a maximum of 7 days after which it must be replaced.

The device can be used only for CPAP (Continuous Positive Airway Pressure) Non Invasive Ventilation treatments.

The **valve functioning at 25 cm H₂O** can be used only with Starmed brand CPAP (Continuous Positive Airway Pressure) Non Invasive Ventilation helmets equipped with a port dedicated to sealed access on the film (Fig. 2).

INSTRUCTIONS FOR USE

DEVICE PREPARATION

Remove the device from the packaging.

Verify the product integrity.

Verify the mobility of the internal plate.

For the use of valve CP918 keeping internal pressure lower than about 20 cm H₂O

Introduce the device into the respiratory system by means of dedicated connections.

For the use of valve keeping internal pressure \leq 25 cm H₂O

Introduce the device into one of the ports as indicated below:


Fig. 3 Sealed access into the front port of Starmed helmet


Fig. 4 Pressure relief valve to be introduced into the front port of Starmed helmet

Remove completely the sealed access from the helmet front port, if present. The front port must without any sealed accesses.

Connect the relief valve into the helmet front port up to the end of the connection (Fig. 4).

The device is ready to be used, as well as to administer CPAP ventilation gases into the helmet.

DEVICE DISCONNECTION

In order to remove the relief valve from the helmet front port, it is necessary to pull on the side tabs and, at the same time, remove the device.

Before administering again CPAP ventilation gases, introduce again the sealed access into the front port.

WARNINGS AND PRECAUTIONS

The activation of the pressure relief valve, indicates that the 20 cmH₂O pressure threshold has been crossed within the breathing circuit.

Intersurgical declines all responsibility for conformity of the product if improperly reused by the customer, given that:

- The physical configuration of the product does not allow adequate cleaning and decontamination with removal of biological residues possibly present on the used device;
- Insufficient data are available to assess the chemical/toxicological safety of the product after exposure of the materials to disinfecting agents;

- Insufficient data are available to assess the maintenance of performance and physical safety of the product after its exposure to cleaning and disinfecting agents;
- Product performance over time has been tested for the life cycle established by the manufacturer.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.

2. Caution: The device contains metal; do not use in departments where the presence of metal is a source of risk for the safety of the patient or third parties.

3. The device can be used for maximum 7 days. Depending on the clinical condition of the patient the medical/nursing staffs is responsible for an eventual more frequent replacement of the device.

4. Non-sterile. Do not sterilize.

5. Do not use on another patient.

6. Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).

7. Dispose of the materials immediately after use according to the applicable regulations.

8. In order to function correctly, the relief valve requires a minimum physiological output of ventilation flow during the treatment. This leakage is not significantly affected by variation of ventilation flows, it does not affect the system functioning and therefore it does not affect the CPAP pressure therapy.

Valve in Fig. 1 allows not to exceed 20 cm H₂O of the helmet internal pressure. It can be used only with devices with connections compliant to standard ISO 5356-1.

Valve in Fig. 2 allows not to exceed 25 H₂O of the helmet internal pressure.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Make reference also to IFU attached to NIV helmets.

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