




## Incentive spirometry

SPIRO-BALL® incentive exerciser and TRI-BALL® respiratory exerciser



Oxygen and Aerosol Therapy ▪ Respiratory Physiotherapy

## SPIRO-BALL®

SPIRO-BALL® Incentive exerciser is a single patient use device designed to assist patients who have undergone abdominal or thoracic intervention and for patients with diseases associated with a decrease in respiratory capacity to recover normal breathing. The device provides a visual aid for both patient and clinician to indicate and monitor progress. The clinician can set a volume target for the patient to try and reach whilst using the device.

### Flow indicator

Intuitive flow control with a smiling face



### Target pointer

Promote the exercise with the target pointer

### Inspiration tube

Flexible tube with mouthpiece

### Two particle filters in the inspiration system

Minimise the risk of inhaling foreign bodies

Code	Description	Box Qty.
1537000	SPIRO-BALL® Incentive Exerciser 4000ml (adult)	12
1539000	SPIRO-BALL® Incentive Exerciser 2500ml (paediatric)	12

## TRI-BALL®

TRI-BALL® encourages the patient to achieve a maximum prolonged inhalation, named SMI technique (Sustained Maximal Inspiration). The SMI is an exercise of deep, slow breathing that when repeated several times per day, may reduce respiratory complications.

### Flexible tube with mouthpiece

Includes a particle filter



### Single patient use

Reduces the risk of cross-contamination

### Three coloured balls in independent chambers

Allow exercise of three different levels of inspiration

Code	Description	Box Qty.
1538000	TRI-BALL® Respiratory Exerciser	12

Both products can be used for up to 30 days.

Please note these products are only available in select markets. For further information, contact your local office.

SPIRO-BALL® and TRI-BALL® are registered trademarks of Plastic Electromechanic Company.

IS10.14 • Issue 8 04.26