



Hydro-Guard™ Pro breathing filter



Airway Management ▪ Breathing filters

Hydro-Guard™ Pro breathing filter

The Hydro-Guard Pro is a low volume pleated mechanical breathing filter which provides high efficiency HEPA filtration. The pleated mechanical filter media provides high efficiency bacterial and viral protection of $\geq 99.9999\%$ and the low volume minimises potential re-breathing, making it ideal for use at the patient airway and is available as a straight or angled option.

Tethered luer cap

For safe monitoring of respiratory gases

Pleated mechanical filter

Provides high efficiency HEPA filtration

Straight and angled options

To meet all clinical needs

Colour coding

Enhances identification, ensuring that the correct filter/HME is used for the correct clinical application



Hydro-Guard Pro

Hydro-Guard Pro Angled

Code	1745162	1745152
Box Qty.	40	40
Filtration efficiency	99.9999%	99.9999%
Resistance at 30 L/min	1.8 cm H ₂ O	1.8 cm H ₂ O
Resistance at 60 L/min	3.8 cm H ₂ O	3.8 cm H ₂ O
Compressible volume	23 ml	25 ml
Weight	23 g	23 g
Connectors	22M/15F-15M	22M/15F-15M
Minimum tidal volume	≥ 150 ml	≥ 150 ml

Make an enquiry

IS12.29 • Issue 3 07.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.