




## Anaesthetic face masks without a hook ring

A clear choice to reduce your impact on the environment




**Anaesthesia** ▪ Anaesthetic Face Masks

## ClearLite™ anaesthetic face masks

A range of anaesthetic face masks which includes seven sizes to meet the clinical needs of all patient types and sizes. The range is available in sizes 0 neonatal to size 6 large adult.

Due to the enhanced anatomical design, the ClearLite mask offers superior seal performance with minimum pressure, reducing the risk of unnecessary discomfort or potential

trauma. The superior seal will also assist the clinician in providing optimum manual ventilation of the patient in both the anaesthetic and emergency situation.

ClearLite is non-PVC and phthalate free which significantly minimises the impact upon the environment. 

### Lightweight mask

Provides enhanced grip and comfort in use for the clinician

### Anatomical shape

Provides a precise fit for different face types and sizes. Minimum pressure required to achieve an effective seal reducing the risk of inadvertent trauma and nerve damage

### Crystal clear body

Provides excellent visibility for confirmation that the patient is well perfused

### Flexible non-inflatable cuff

Provides a comfortable and effective airtight seal



### Safe and secure connection

This is achieved between the mask and the breathing system in use. The ISO 22 mm tapered connection complies with ISO 5356

### Reduced mask volume

Reduces potential rebreathing and risks associated with elevated levels of end tidal CO<sub>2</sub>

### Colour-coding

The comfortable, flexible seal is colour coded to allow quick and easy identification of the correct sized mask

### Low environmental impact

Non-PVC and phthalate free



Code	Description	Colour	Box Qty.
7296001	ClearLite, anaesthetic face mask, size 6, extra large adult, red seal, 22F		20 
7295001	ClearLite, anaesthetic face mask, size 5, large adult, orange seal, 22F		35 
7294001	ClearLite, anaesthetic face mask, size 4, adult, green seal, 22F		35 
7293001	ClearLite, anaesthetic face mask, size 3, small adult, yellow seal, 22F		35 
7292001	ClearLite, anaesthetic face mask, size 2, paediatric, white seal, 22F		25 
7291000	ClearLite, anaesthetic face mask, size 1, infant, grey seal, no hook ring, 15M		30 
7290000	ClearLite, anaesthetic face mask, size 0, neonate, blue seal, no hook ring, 15M		40 

[Make an enquiry](#)

[Find out more](#)

[Watch the video](#)



## Economy anaesthetic face masks

A range of seven colour-coded economy PVC anaesthetic face masks with a soft cushion and clear housing. Seven sizes from neonatal to extra large adult.



22F

### Clear flexible housing and cushion

For good patient visibility

### Soft cushion

For patient comfort and air tight fit

### Seven sizes available

From neonatal to extra large adult



Code	Description	Colour	Box Qty.
1517111	Economy, anaesthetic face mask, size 6, extra large adult, 22F	Red	25
1516111	Economy, anaesthetic face mask, size 5, large adult, 22F	Orange	30
1515111	Economy, anaesthetic face mask, size 4, medium adult, 22F	Green	40
1514111	Economy, anaesthetic face mask, size 3, small adult, 22F	Yellow	50
1513111	Economy, anaesthetic face mask, size 2, paediatric, 22F	White	35
1512111	Economy, anaesthetic face mask, size 1, infant, 15M	Black	40
1511111	Economy, anaesthetic face mask, size 0, neonate, 15M	Blue	50

Make an enquiry

## Scented anaesthetic face masks

### In a range of sizes and scents

A range of scented single use PVC anaesthetic face masks with adjustable cushion, available in seven sizes and three scents: cherry, strawberry or vanilla, presented in a cool mint colour to help the patient relax.



22F

#### Cushion valve at top (nose) of face mask

For user comfort

#### Soft adjustable face cushion

For patient comfort and a good anatomical fit



**Clear flexible housing**  
For good patient facial visibility

**Scented**  
To help patient relax

Code	Description	Colour	Box Qty.
1126111	Scented, vanilla, anaesthetic face mask, size 6, extra large adult, 22F	●	24
1125111	Scented, vanilla, anaesthetic face mask, size 5, large adult, 22F	●	30
1124111	Scented, vanilla, anaesthetic face mask, size 4, medium adult, 22F	●	40
1123111	Scented, vanilla, anaesthetic face mask, size 3, small adult, 22F	●	50
1122111	Scented, cherry, anaesthetic face mask, size 2, paediatric, 22F	○	35
1129111	Scented, strawberry, anaesthetic face mask, size 2, paediatric, 22F	○	35
1121111	Scented, cherry, anaesthetic face mask, size 1, infant, 15M	●	32
1128111	Scented, strawberry, anaesthetic face mask, size 1, infant, 15M	●	32
1120111	Scented, cherry, anaesthetic face mask, size 0, neonate, 15M	●	40
1127111	Scented, strawberry, anaesthetic face mask, size 0, neonate, 15M	●	40

Make an enquiry

IS5.35 • Issue 2 09.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.