



## OptiLube™ sachets of lubricant

For use with supraglottic airways, endotracheal tubes and naso/oropharyngeal airways



**Airway Management** ▪ Airway Accessories

## OptiLube™ sachets of lubricant

### For use with supraglottic airways, endotracheal tubes and naso/oropharyngeal airways

To complement our wide range of airway devices, we are delighted to announce the addition of sachets of lubricant to our product portfolio as a service to our customers.

OptiLube™ sachets of lubricant contain 5g of water soluble lubricating jelly for use in the preparation of airway devices such as supraglottic airways, endotracheal tubes and naso/oropharyngeal airways. It can also be used to prime the gastric channel of second generation airway devices such as the i-gel® or lubricate suction and orogastric tubes.

OptiLube sachets are harmless to human tissue and will not injure rubber, metal appliances, instruments and natural or synthetic textiles.



Code	Description	Box Qty.
5469000	OptiLube™, sachets of lubricant, 5g	150

[Make an enquiry](#)

OptiLube is a trademark of Optimum Medical Solutions.



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
 T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.

IS6.8 ▪ Issue 4 07.19